




Creating
Perfectly Sculptured
Soft Tissue...

...for individual esthetics



ProHeal Cap 20, 4.3 mm, 5.5 mm Short and 5.5 mm Long

The ProHeal Cap system allows for sculpturing of the soft tissue to create individual esthetics when making screw-retained prosthetics.

The ProHeal Cap, available in three different sizes offers flexibility and an easy clinical procedure to sculpt the peri-implant soft tissue. The corresponding Abutment Impression Pick-up allows for perfect registration of the soft tissue at impression taking.

Ideal soft tissue sculpturing

The ProHeal Cap is a system of sterile healing caps providing ultimate soft tissue esthetics for screw-retained prosthetics. The system consists of a tissue sculpturing healing cap with

a corresponding Impression Pick-up placed on the Uni Abutment during the healing period. The ProHeal Cap controls and shapes the surrounding soft tissue.

The ProHeal Cap also protects the conical top surface and the screw entrance on the Uni Abutment from plaque and tissue overgrowth. The material, design and surface texture of the healing caps provide smooth feeling, tissue health and patient comfort.

Guide Pin Hex
13 and 16 mm



Guide Pin Hex
19 and 22 mm



Creating perfect impression

The Abutment Impression Pick-up corresponding to the ProHeal Cap simplifies the impression taking and is the key for accurate impressions.

The novel healing and tissue sculpturing system also include Guide Pin Hex.

- The Abutment Impression Pick-up is a two-component pick-up shaped specifically to mimic the ProHeal Cap and utilized for taking impressions on the Uni Abutment.
- The Guide Pin Hex represents a new type of guide pin with a knurled top and hexagonal head for improved handling.

For further information on the ProHeal Cap System, contact your local Astra Tech representative.



The ProHeal Caps in place.



Sculptured papillas and soft tissue after healing.



The corresponding Abutment Impression Pick-ups installed.

The Abutment
Impression
Pick-up 20
4.3 mm



The Abutment
Impression
Pick-up 20
5.5 mm



ASTRA TECH IMPLANTS

D E N T A L S Y S T E M

The Astra Tech Dental Implant System is based on years of scientific research and clinical documentation. The result is a unique system with perfect balance of **Simplicity, Reliability** and **Esthetics**, proven in long-term clinical studies. Astra Tech is committed to forming long-term working partnerships with dental professionals, offering products and services that improve quality of life.

Denmark

Astra Tech A/S, Husby Allé 19, DK-2630 Taastrup.
Tel: +45 43 71 33 77. Fax: +45 43 71 78 65.

Finland

Astra Tech Oy, PL 96, FI-02230 Espoo.
Tel: +358 9 8676 1626. Fax: +358 9 8044 128.

France

Société de Matériel Médical, Astra Tech France,
Groupe Pharmaceutique AstraZeneca, 1, rue de l'Union,
TSA 90002, FR-92843 Rueil-Malmaison Cedex.
Tel: +33 1 41 39 02 40. Fax: +33 1 41 39 02 44.

Germany

Astra Tech GmbH, An der kleinen Seite 8, DE-65604 Elz.
Tel: +49 6431 9869 0. Fax: +49 6431 9869 30.

Italy

Astra Tech S.p.A
Via M. Luther King n. 38/2, IT-40132 Bologna
Tel: +39 051 641 6511. Fax: +39 051 641 7742.

The Netherlands

Astra Tech BV, Signaalrood 55, NL-2718 SG Zoetermeer.
Tel: +31 79 360 1950. Fax: +31 79 362 3748.

Norway

Astra Tech AS, Postboks 160, NO-1471 Lørenskog.
Tel: +47 67 92 05 50. Fax: +47 67 92 05 60.

Spain

Astra Tech S.A.,
Calle Ciencias nº 73 derecha. Nave 9, Polígono
Industrial Pedrosa, ES-08908 L'Hospitalet de Llobregat.
Tel. Servicio al cliente: +34.902.101.558.
Tel: +34.932.643.560. Fax: +34.933.363.231.

Sweden

Astra Tech AB, P.O. Box 14, SE-431 21 Mölndal.
Tel: +46 31 776 30 00. Fax: +46 31 776 30 10.

United Kingdom

Astra Tech Ltd., Brunel Way, Stonehouse, Glos GL10 3SX.
Tel: +44 1453 791763. Fax: +44 1453 791001.

USA

Astra Tech Inc.
430 Bedford Street, Suite 100, Lexington, MA 02420.
Tel: +1 781 861 7707. Fax: +1 781 861 7787.

Japan

Astra Tech Division, AstraZeneca K.K.
Tokyo Regional Office
Koraku Mori Bldg. 18F, 1-4-14, Koraku
Bunkyo-ku, Tokyo 112-0004.
Tel: +81 3 5840 1113. Fax: +81 3 5840 1151.

Other Markets

Astra Tech AB, Export Department,
P.O.Box 14, SE-431 21 Mölndal, Sweden.
Tel: +46 31 776 30 00. Fax: +46 31 776 30 23.

